

## Do's and Don'ts

### Do's

- This product should always place on plane surface.
- Always plug in or plug out Power connector with dry hands.
- Use caution when raising or folding the stand of the equipment to avoid injury. When the stand is raised, do not touch the joints or any other moving part of the stand. Moving parts.

### Don'ts

- Don't open BMI stand in open place when air is blowing fast.
- To prevent electrical shock or damage device don't install BMI machine on wet place and connect power adopter to appropriate power supply (230v AC +/- 10% and 50 Hz +/- 2%).
- Don't jump on weight machine and don't hold stand during measuring in order to measure accurate result.
- Don't open the device. Internal parts are not for customer use. If the equipment is open, the warranty is void, and service costs will be charged. If service is required, contact to our company or the supplying agency.
- Don't touch the sensor and don't change sensor cover.



### WARNING

Physically disabled person or children should avoid measure BMI to avoid injury or measure under assistance.